COULTER 5C Cell Control

TABLE OF CONTENTS

| | |
|-----------------------|----|
| English | 2 |
| Français | 5 |
| Deutsch | 8 |
| Italiano | 11 |
| Español | 14 |
| Português | 17 |
| Dansk | 20 |
| Svenska | 23 |
| Norsk | 26 |
| Ελληνικά | 29 |
| 日本語 | 32 |
| Lietuviškai | 35 |
| Magyar | 38 |
| Polski | 41 |
| Čeština | 44 |
| 한국어 | 47 |
| Türkçe | 50 |
| Русский | 53 |
| Hrvatski | 56 |

COULTER 5C Cell Control

REF 7547199, 7547124, 7547117,
7547116, 7547011, 7547010,
7547001, 624506

**For *In Vitro* Diagnostic Use
Rx Only**

INTENDED USE

5C Cell Control is a hematology quality control material used to monitor the performance of COULTER hematology analyzers listed in the TABLE OF EXPECTED RESULTS in conjunction with specific COULTER reagents. Refer to your Product Manuals or Online Help System.

The assigned values and expected ranges on the TABLE OF EXPECTED RESULTS can be used to monitor instrument performance. This product can also be used to establish your own laboratory mean.

SUMMARY AND PRINCIPLE

5C Cell Control is a reference product prepared from stabilized human blood. By design, 5C Cell Control confirms and monitors instrument accuracy and precision performance by providing measurements for counting, sizing, hemoglobin determination and White Blood Cell differentiation using VCS technology.

REAGENTS

5C Cell Control consists of treated, stabilized human erythrocytes in an isotonic medium. 5C Cell Control also contains a stabilized, platelet-sized component, and fixed erythrocytes to simulate leukocytes.

WARNING AND PRECAUTIONS

POTENTIAL BIOHAZARDOUS MATERIAL

Product contains biologically sourced materials: Human, Avian, Reptile and Ungulate.

Each human donor unit used in preparation of this material was tested by an FDA approved method for the presence of the antibodies to Human Immunodeficiency Virus (HIV-1 and HIV-2) and Hepatitis C Virus (HCV) as well as for hepatitis B virus surface antigen and found to be negative (were not repeatedly reactive).

Because no test method can offer complete assurance that hepatitis B virus, Human Immunodeficiency Virus (HIV-1 and HIV-2), or other infectious agents are absent, this specimen/reagent should be handled at Biosafety Level 2, as recommended for any potentially infectious human serum or blood specimen in the Centers for Disease Control/National Institutes of Health manual "Biosafety in Microbiological and Biomedical Laboratories," 1988.

Product contains <0.1% Sodium Azide. Sodium Azide preservative may form explosive compounds in metal drain lines. See National Institute for Occupational Safety and Health Bulletin: Explosive Azide Hazards (8/16/76).

| GHS HAZARD CLASSIFICATION | | |
|---|--|-----------------------------------|
| 5C CELL CONTROL, Abnormal I, Abnormal II and Normal | | May produce an allergic reaction. |
| | | Neomycin Sulfate <0.1% |
| | | Penicillin-G, sodium salt <0.1% |

STORAGE, STABILITY, AND DISPOSAL

5C Cell Control is shipped in a thermally insulated container designed to keep it cool. When stored at 2 to 8°C, sealed/unopened tubes are stable until the expiration date shown on the TABLE OF EXPECTED RESULTS. For opened vial stability, refer to the TABLE OF EXPECTED RESULTS for your system.

NOTE: The MCV, RDW and/or RDW-SD parameters may show trending through the product's shelf life. This is inherent to the product and should not be considered an indicator of product instability. 95% of recovered values for these parameters should remain within the stated EXPECTED RANGES.

Storage of product with cap down (inverted) may require additional mixing for complete resuspension of cellular components.

Dispose of waste product, unused product and contaminated packaging in compliance with federal, state and local regulations.

EVIDENCE OF DETERIORATION

Inability to obtain expected values in the absence of known instrument problems or gross hemolysis (darkly colored supernatant) is indicative of product deterioration. However, a slight pink color to the supernatant is normal and should not be confused with deterioration of the product.

MATERIALS PROVIDED

- 5C Cell Control. (Refer to the Product Availability section of this insert for quantities and fill volumes for your particular kit).
- Table of Expected Results containing Assigned Values and Expected Ranges.
- Removable Media, if applicable.

MATERIALS REQUIRED, BUT NOT PROVIDED

- Removable Media, if applicable.
- Sufficient COULTER reagents to complete analysis.
- Absorbent lint-free material such as paper wipes.
- Product Manuals or Online Help System.

INSTRUCTIONS FOR USE

1. Remove 5C Cell Control tube(s) from refrigerator and warm at AMBIENT TEMPERATURE for 10 to 15 minutes.
2. Prepare your instrument: For COULTER LH 500, LH 750, LH780 and HmX systems: Insert Removable Media with Assay Values when a new lot of 5C Cell Controls is received. Follow the instructions on the Workstation/computer screen, and in the Online Help or Product Manuals.
3. After warming, mix by hand as follows:
 - A. Roll the tube slowly between the palms of the hands eight times in an upright position.
 - B. Invert the tube and slowly roll it between the palms eight times.
 - C. Gently invert the tube eight times.

4. Repeat steps 3A to 3C.

NOTE: If your COULTER instrument includes an autosampler/mixer, you must mix first as directed above, then place control tubes into a cassette and place on the instrument. DO NOT use any other type of mechanical mixer.

5. Run the control in the Primary mode of the instrument.

NOTE: For an instrument without cap-piercing capability, analyze the control as instructed in the Product Manuals or Online Help System. For your protection, use an absorbent lint-free material such as paper wipes when removing the cap/stopper assembly from the tube. After use, immediately wipe the cap and mouth of the tube before replacing cap.

6. Return the tube(s) to the refrigerator within **30 minutes**.
7. Compare instrument values to those given in the TABLE OF EXPECTED RESULTS.

ASSIGNED VALUES AND EXPECTED RESULTS

Refer to TABLE OF EXPECTED RESULTS.

For LH 750 and LH 780, the WBC assigned value and expected range provided in the TABLE OF EXPECTED RESULTS only applies to the WBC value; it does not apply to uncorrected WBC (UWBC) value.

ASSIGNED VALUES are determined on validated systems using specific COULTER reagents. The values are confirmed by multiple analysis of the control.

EXPECTED RANGES include variation between lots and between instruments, and represent 95% confidence limits for well maintained instrument systems.

Whether you use the Beckman Coulter Assigned Value or your own laboratory mean, the instrument is considered well maintained and operating correctly if:

- $\geq 95\%$ of the recovered values fall within the EXPECTED RANGE of the Beckman Coulter Assigned Value.
- Recovered values do not trend OUTSIDE the EXPECTED RANGE.

Before your current cell control lot(s) expire, perform the following on your new lot(s):

- Confirm that recovered values are within the TABLE OF EXPECTED RESULTS.
- OR
- Establish your own laboratory mean.

For LH 500, LH 750 and LH780 Systems:

Reviewed flagged ('R') differential results can be accepted for instrument systems where the following conditions have been verified by the operator:

1. Patient results are exhibiting an acceptable rate of Suspect flagging for "Verify Diff".
2. Your instrument is operating in the proper operating temperature mode.
3. 95% of the recovered Control values for the Differential fall within the Expected Range as listed on the Table of Expected Results for the lot number in use, and recovered values do not trend outside the Expected Range.
4. The system is free of System Alarms when 'R' flags are observed with the Control Differential results.

PERFORMANCE LIMITS

Individual laboratories can expect better precision than that shown in the EXPECTED RANGE column. Refer to your Product Manuals or Online Help System for precision performance characteristics for your instrument.

TROUBLESHOOTING PROCEDURE

The Product Manuals or Online Help System contains troubleshooting procedures for resolving control recovery problems. If you need additional help, call Beckman Coulter Customer Service at 800-526-7694 (USA and Canada) or contact your local Beckman Coulter Representative.

To provide faster handling of your inquiry, please have the following information available when you call:

- Instrument model and serial number.
- Expiration dates, lot numbers, and supporting data for all reagents, for control(s) in question and for other Cell Control lots currently in use.
- Previous Cell Control lot numbers and the data for these previous lots.
- Data history for COULTER LATRON control analyses.
- Data from current reproducibility study (N=10) using a fresh whole blood specimen and performed according to the Product Manuals or Online Help System.
- Data from the IQAP survey as well as data from your last instrument calibration, where appropriate.

PRODUCT AVAILABILITY

COULTER 5C Cell Control

| | |
|------------|--|
| REF | 624506 - 1 x 3.3 mL Normal, 1 x 3.3 mL Abnormal I, 1 x 3.3 mL Abnormal II |
| REF | 7547001 - 3 x 3.3 mL Normal, 3 x 3.3 mL Abnormal I, 3 x 3.3 mL Abnormal II |
| REF | 7547010 - 9 x 3.3 mL Abnormal I |
| REF | 7547011 - 9 x 3.3 mL Abnormal II |
| REF | 7547116 - 4 x 3.3 mL Normal, 4 x 3.3 mL Abnormal I, 4 x 3.3 mL Abnormal II |
| REF | 7547117 - 9 x 3.3 mL Normal |
| REF | 7547124 - 4 x 3.3 mL Normal, 4 x 3.3 mL Abnormal I, 4 x 3.3 mL Abnormal II |
| REF | 7547199 - 4 x 3.3 mL Normal, 4 x 3.3 mL Abnormal I, 4 x 3.3 mL Abnormal II |

TRADEMARKS

Beckman Coulter, the stylized logo, COULTER and 5C are trademarks of Beckman Coulter, Inc. and are registered in the USPTO.

For additional information or if damaged product is received, call Beckman Coulter Customer Service at 800-526-7694 (USA or Canada) or contact your local Beckman Coulter Representative.

ADDITIONAL INFORMATION

SDS Safety Data Sheet is available at techdocs.beckmancoulter.com

Revision History

Revision FE, 12/2012

Changes were made to:

- Added Korean language

Revision FF, 10/2014

Changes were made to:

- GHS HAZARD CLASSIFICATION
- Added Croatian Language

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Beckman Coulter Ireland Inc.
Beckman Coulter, Inc., 250 S. Kraemer Blvd., Brea, CA 92821 U.S.A.